

Pancreatic islet culture medium Kit

Catalog No: M2269-Kit includes:

- 5.0 ML L-Glutamine
- 5.0 ML Penicillin-Streptomycin Solution
- 25.0 ML FBS

*Before use, add M2269-Kit to 500ml basal medium.

Storage Condition

Store the basal medium at 2-8°C. Store cell growth supplement, fetal bovine serum (FBS) and antibiotics at -20°C. The complete cell culture medium with Supplement Kit can be kept in 4°C for two months. Protect from light.

To assure sterility, the prepared medium should be filtered with a 0.2 um filter after 2 weeks or if there is concern that sterility was compromised during the supplementation process.

Shipping

Ice pack and dry ice.

Authorized Uses of *Cell Biologics* Products

Pancreatic islet culture medium Kit from *Cell Biologics* are distributed for internal research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' cells or products from the purchaser to other markets, organizations or individuals is prohibited by *Cell Biologics*, without the company's written consent. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.