

Leukocyte Medium Kit

Catalog Number: M3366-Kit

Description

Leukocyte Medium Supplement Kit, Cat. No. M3366-Kit includes:

- 5.0 ML Penicillin-Streptomycin
- 50.0 ML FBS

*Before use, add M3366-Kit to 500ml basal medium.

Storage Condition

Store endothelial cell growth supplement, fetal bovine serum (FBS) and antibiotics at -20°C. The complete cell culture medium with Supplement Kit can be kept in 4°C for two months. Protect from light.

Note: To assure sterility after 2 weeks or if there is concern that sterility was compromised during the supplementation process, the prepared medium may be refiltered with a 0.2 um filter.

Shipping

Dry Ice.

Authorized Uses of *Cell Biologics* Products

Leukocyte Medium Kit from *Cell Biologics* is distributed for internal research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' Cells or Products from the purchaser to other markets, organizations, or individuals is prohibited by *Cell Biologics*. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.