

## Neutrophil Medium Kit

Catalog No: M3367-Kit

### Description

Neutrophil Medium Supplement Kit, Cat. No. M3367-Kit includes:

- 5.0 ML MEM Non-Essential Amino Acids Solution (100x)
- 5.0 ML Penicillin-Streptomycin
- 50.0 ML FBS

\*Before use, add M3367-Kit to 500ml basal medium.

### Storage Condition

Store cell growth supplement, fetal bovine serum (FBS) and antibiotics at -20°C. The complete cell culture medium with Supplement Kit can be kept in 4°C for two months. Protect from light.

To assure sterility, the prepared medium should be filtered with a 0.2 um filter after 2 weeks or if there is concern that sterility was compromised during the supplementation process.

### Shipping

Dry Ice.

### Authorized Uses of *Cell Biologics* Products

Neutrophil Medium Kit from *Cell Biologics* is distributed for internal research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' cells or products from the purchaser to other markets, organizations or individuals is prohibited by *Cell Biologics* without the company's written consent. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.