

**Catalog Number: M6621-Kit**

**Description**

Epithelial Cell Medium Supplement Kit, Cat. No. M6621-Kit includes:

- 0.5 ML Hydrocortisone
- 0.5 ML EGF
- 5.0 ML Antibiotic-Antimycotic Solution
- 10.0 ML FBS

\*Before use, add M6621 Supplement Kit (M6621-Kit) into 500ml basal medium.

**Storage Condition**

Store endothelial cell growth supplement, fetal bovine serum (FBS) and antibiotics at -20°C. The complete cell culture medium with Supplement Kit can be kept in 4°C for two months. Protect from light.

Note: To assure sterility after 2 weeks or if there is concern that sterility was compromised during the supplementation process, the prepared medium may be refiltered with a 0.2 um filter.

**Shipping**

Dry Ice.

**Authorized Uses of Cell Biologics Products**

Epithelial Cell Medium Kit from *Cell Biologics* is distributed for in vitro research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' Cells or Products from the purchaser to other markets, organizations, or individuals is prohibited by Cell Biologics. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.