

Human Endothelial Cell Medium Kit

Catalog Number: H1168-Kit

Description

Endothelial Cell Medium Supplement Kit, Cat. No. H1168-Kit includes:

- 0.5 ML VEGF
- 0.5 ML Heparin
- 0.5 ML EGF
- 0.5 ML FGF
- 0.5 ML Hydrocortisone
- 5.0 ML L-Glutamine
- 5.0 ML Antibiotic-Antimycotic Solution
- 10.0 ML Endothelial Cell Supplement
- 50.0 ML FBS

*Before use, add H1168 Supplement Kit (H1168-Kit) into 500ml basal medium.

Storage Condition

Store the basal medium at 2-8°C. Store endothelial cell growth supplement, fetal bovine serum (FBS) and antibiotics at -20°C. The complete cell culture medium with Supplement Kit can be kept in 4°C for two months. Protect from light.

Note: To assure sterility after 2 weeks or if there is concern that sterility was compromised during the supplementation process, the prepared medium may be refiltered with a 0.2 um filter.

Shipping

Ice pack and dry ice.

Authorized Uses of *Cell Biologics* Products

Human Endothelial Cell Medium Kit from *Cell Biologics* are distributed for internal research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' Cells or Products from the purchaser to other markets, organizations, or individuals is prohibited by *Cell Biologics*. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.