

## Endothelial Cell Medium Kit

Catalog Number: M1168-Kit

### Description

Endothelial Cell Medium Supplement Kit, Cat. No. M1168-Kit includes:

- 0.5 ML VEGF
- 0.5 ML ECGS
- 0.5 ML Heparin
- 0.5 ML EGF
- 0.5 ML Hydrocortisone
- 5.0 ML L-Glutamine
- 5.0 ML Antibiotic-Antimycotic Solution
- 25.0 ML FBS

\*Before use, add M1168 Supplement Kit (M1168-Kit) into 500ml basal medium.

### Storage Condition

Store the basal medium at 2-8°C. Store endothelial cell growth supplement, fetal bovine serum (FBS) and antibiotics at -20°C. The complete cell culture medium with Supplement Kit can be kept in 4°C for two months. Protect from light.

Note: To assure sterility after 2 weeks or if there is concern that sterility was compromised during the supplementation process, the prepared medium may be refiltered with a 0.2 um filter.

### Shipping

Ice pack and dry ice.

### Authorized Uses of *Cell Biologics* Products

Endothelial Cell Medium Kit from *Cell Biologics* are distributed for internal research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' Cells or Products from the purchaser to other markets, organizations, or individuals is prohibited by *Cell Biologics*. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.