

## Complete Fibroblast Medium /w Kit – 500 ML

Catalog Number: M2267  
(500 ml Basal medium with the growth factor supplement)

### Description

Fibroblast Medium is a complete medium designed for the culture of fibroblasts. It was tested and optimized with fibroblast growth and proliferation *in vitro*. It is formulated for use with 5% CO<sub>2</sub> and 95% air in a humidified incubator.

The medium consists of 500 ml of basal medium (containing essential and non-essential amino acids, vitamins, organic and inorganic compounds, hormones, growth factors, trace minerals), supplemented with cell growth supplement, antibiotics, and fetal bovine serum.

Fibroblast Medium Supplement Kit, Cat. No. M2267-Kit includes:

- 0.5 ML FGF
- 0.5 ML Hydrocortisone
- 5.0 ML L-Glutamine
- 5.0 ML Antibiotic-Antimycotic Solution
- 50.0 ML FBS

\*Before use, add M2267 Supplement Kit (M2267-Kit) into 500ml basal medium.

### Storage Condition

Store the basal medium at 2-8°C. Store cell growth supplement, fetal bovine serum (FBS) and antibiotics at -20°C. The complete cell culture medium with Supplement Kit can be kept in 4°C for two months. Protect from light.

Note: To assure sterility after 2 weeks or if there is concern that sterility was compromised during the supplementation process, the prepared medium may be refiltered with a 0.2 um filter.

### Shipping

Ice pack and dry ice.

### Authorized Uses of *Cell Biologics* Products

Fibroblast Media from *Cell Biologics* are distributed for *in vitro* research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' Cells or Products from the purchaser to other markets, organizations, or individuals is prohibited by *Cell Biologics*. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.