

**Complete Macrophage Medium /w Kit – 500 ML**

Catalog No: M3368  
(500 ml Basal medium with M3368-Kit)

**Description**

Macrophage Medium is a complete medium designed for the culture of human and animal macrophages. It was tested and optimized with human and animal macrophage growth *in vitro*. It is formulated for use with 5% CO<sub>2</sub> and 95% air in a humidified incubator.

The medium consists of 500 ml of basal medium (containing essential and non-essential amino acids, vitamins, organic and inorganic compounds, hormones, growth factors, trace minerals), supplemented with murine macrophage growth supplement, antibiotics, and fetal bovine serum.

Macrophage Medium Supplement Kit, Cat. No. M3368-Kit includes:

- 40 ML GM-CSF
- 50 ML FBS
- 5.0 ML Penicillin-Streptomycin

\*Before use, add M3368-Kit to 500ml basal medium.

**Storage Condition**

Store the basal medium at 2-8°C. Store cell growth supplement, fetal bovine serum (FBS) and antibiotics at -20°C. The complete cell culture medium with Supplement Kit can be kept in 4°C for two months. Protect from light.

To assure sterility, the prepared medium should be filtered with a 0.2 um filter after 2 weeks or if there is concern that sterility was compromised during the supplementation process.

**Shipping**

Ice pack and dry ice.

**Authorized Uses of *Cell Biologics* Products**

Macrophage Media from *Cell Biologics* are distributed for internal research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' cells or products from the purchaser to other markets, organizations or individuals is prohibited by *Cell Biologics*, without the company's written consent. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.