

**Macrophage Culture Supplement – 40.0 ML**

Catalog Number: CB3368

**Description**

Macrophage Culture Supplement (L-cell supernatant) is used to cultivate macrophages. Usually it can be used in medium for 7 days.

**Storage Condition**

Store the solution at -20°C.

**Shipping**

Dry Ice.

**Authorized Uses of *Cell Biologics* Products**

Macrophage Culture Supplements from *Cell Biologics* are distributed for internal *in vitro* research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' cells or products from the purchaser to other markets, organizations or individuals is prohibited by *Cell Biologics*, without the company's written consent. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.