

Gelatin-Based Coating Solution

Catalog Number: CB6950

Description

Gelatin-Based Coating Solution is designed to coat cell culture flasks or dishes. Coat sterile culture dishes or flasks with Gelatin-Based Coating Solution for 2 min and then aspirate the excess solution before seeding cells. The use of coated-flasks improves cell attachment for certain types of primary cells.

Storage Condition

Store the solution at 2-8°C.

Shipping

Ice pack.

Authorized Uses of *Cell Biologics* Products

Gelatin-Based Coating Solution from *Cell Biologics* are distributed for internal *in vitro* research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' cells or products from the purchaser to other markets, organizations or individuals is prohibited by *Cell Biologics*, without the company's written consent. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.