

CB Cell Lysis Buffer (1X)

Catalog Number: CBL6822

Description

A convenient, ready-to-use cell lysis buffer that is generally used for western blot, Co-IP, Chromatin IP (ChIP), enzyme activity assays, etc.

Lysis buffer contains 50 mM Hepes, 50 mM NaCl, 1% Triton X-100, 5 mM EDTA, 1.5 mM DTT, 10 mM sodium pyrophosphate, 50 mM sodium fluoride, 1 mM sodium orthovanadate, and freshly add protease inhibitor cocktail.

The buffer is sterilized and ready to use for all types of human or animal cells.

Storage Condition

Store the solution at 2-8°C.

Shipping

Ice pack.

Authorized Uses of *Cell Biologics* Products

CB Cell Lysis Buffer (1X) from *Cell Biologics* are distributed for internal *in vitro* research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' cells or products from the purchaser to other markets, organizations or individuals is prohibited by *Cell Biologics*, without the company's written consent. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.