

---

## Serum-Free Endothelial Cell Medium /w Kit – 500 ML

Catalog Number: M1168SF  
(500 ml Basal medium with the growth factor supplement)

### Description

Serum-Free Endothelial Cell Medium is a complete medium designed for the culture of endothelial cells. It was tested and optimized with endothelial cell growth and proliferation *in vitro*. It is formulated for use with 5% CO<sub>2</sub> and 95% air in a humidified incubator.

The medium consists of 500 ml of basal medium (containing essential and non-essential amino acids, vitamins, organic and inorganic compounds, hormones, growth factors, trace minerals), supplemented with endothelial cell growth supplement, and antibiotics.

Endothelial Cell Medium Serum-Free Supplement Kit, Cat. No. M1168SF-Kit includes:

- 0.5 ML VEGF
- 0.5 ML Heparin
- 0.5 ML ECGS
- 0.5 ML EGF
- 0.5 ML Hydrocortisone
- 5.0 ML Antibiotic-Antimycotic Solution

\*Before use, add M1168SF Supplement Kit (M1168SF-Kit) into a 500ml basal medium.

### Storage Condition

Store the basal medium at 2-8°C. Store endothelial cell growth supplement, and antibiotics at -20°C. The complete cell culture medium with Supplement Kit can be kept in 4°C for two months. Protect from light.

Note: To assure sterility after 2 weeks or if there is concern that sterility was compromised during the supplementation process, the prepared medium may be refiltered with a 0.2 um filter.

### Shipping

Ice pack and dry ice.

### Authorized Uses of Cell Biologics Products

Endothelial Cell Media Serum-Free from *Cell Biologics* are distributed for internal research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' Cells or Products from the purchaser to other markets, organizations, or individuals is prohibited by *Cell Biologics*. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.