
Complete Skeletal Muscle Cell Medium (No Phenol Red) /w Kit – 500 ML

Catalog No: M1261PF
(500 ml Basal medium with M1261PF-Kit)

Description

Skeletal Muscle Cell Medium (No Phenol Red) is a complete medium designed for the culture of Skeletal Muscle Cells. It was tested and optimized with Skeletal Muscle Cell growth and proliferation *in vitro*. It is formulated for use with 5% CO₂ and 95% air in a humidified incubator.

The medium consists of 500 ml of basal medium (containing essential and non-essential amino acids, vitamins, organic and inorganic compounds, hormones, growth factors, trace minerals), supplemented with cell growth supplement, antibiotics, and fetal bovine serum.

Skeletal Muscle Cell Medium (No Phenol Red) Supplement Kit, Cat. No. M1261PF-Kit includes:

- 0.5 ML Insulin
- 0.5 ML EGF
- 0.5 ML FGF
- 0.5 ML Hydrocortisone
- 5.0 ML Antibiotic-Antimycotic Solution
- 50.0 ML FBS

*Before use, add M1261PF-Kit to 500ml basal medium.

Storage Condition

Store the basal medium at 2-8°C. Store cell growth supplement, fetal bovine serum (FBS) and antibiotics at -20°C. The complete cell culture medium with Supplement Kit can be kept in 4°C for two months. Protect from light.

To assure sterility, the prepared medium should be filtered with a 0.2 um filter after 2 weeks or if there is concern that sterility was compromised during the supplementation process.

Shipping

Ice pack and dry ice.

Authorized Uses of *Cell Biologics* Products

Skeletal Muscle Cell Medium (No Phenol Red) from *Cell Biologics* are distributed for internal research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' cells or products from the purchaser to other markets, organizations or individuals is prohibited by *Cell Biologics*, without the company's written consent. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.