

Skeletal Muscle Cell Basal Medium – 500 ML

Catalog Number: M1261b

Description

Skeletal Muscle Cell Basal Medium is designed for the culture of human and animal skeletal muscle cells. It was tested and optimized with cell growth and proliferation *in vitro*. It is formulated for use with 5% CO₂ and 95% air in a humidified incubator. The basal medium contains essential and non-essential amino acids, vitamins, organic and inorganic compounds, hormones, growth factors, trace minerals.

Storage Condition

Store the basal medium at 2-8°C for 8 months. Protect from light.

Shipping

Ice pack.

Authorized Uses of *Cell Biologics* Products

Cell Culture Medium from *Cell Biologics* is distributed for internal research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' Cells or Products from the purchaser to other markets, organizations, or individuals is prohibited by *Cell Biologics*. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.