

---

## Complete Cell Culture Medium /w Kit – 500 ML

Catalog No: M5567 (500 ml Basal medium with M5567-Kit)

### Description

Complete Cell Culture Medium is a complete medium designed for the culture of animal and human Cells. It was tested and optimized with Cell culture in vitro. It is formulated for use with 5% CO<sub>2</sub> and 95% air in a humidified incubator.

The medium consists of 500 ml of basal medium (containing essential and non-essential amino acids, vitamins, organic and inorganic compounds, hormones, growth factors, trace minerals), supplemented with cell culture supplement, antibiotics, and fetal bovine serum.

Complete Cell Medium Supplement Kit, Cat. No. M5567-Kit includes:

- 0.5 ML EGF
- 0.5 ML Insulin
- 5.0 ML Hydrocortisone
- 5.0 ML Penicillin-Streptomycin Solution
- 25.0 ML FBS

\*Before use, add M5567-Kit to 500ml basal medium.

### Storage Condition

Store the basal medium at 2-8°C. Store cell growth supplement, fetal bovine serum (FBS) and antibiotics at -20°C. The complete cell culture medium with Supplement Kit can be kept in 4°C for two months. Protect from light.

To assure sterility, the prepared medium should be filtered with a 0.2 um filter after 2 weeks or if there is concern that sterility was compromised during the supplementation process.

### Shipping

Ice pack and dry ice.

### Authorized Uses of *Cell Biologics* Products

Complete Cell Medium from *Cell Biologics* is distributed for internal *in vitro* research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' cells or products from the purchaser to other markets, organizations or individuals is prohibited by *Cell Biologics*, without the company's written consent. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.