

Complete Macrophage Medium /w Kit – 125 ML

Catalog No: M6030

Description

Macrophage Medium is a complete medium designed for the culture of human and animal (Monkey, Porcine, Canine, Rabbit) macrophages. It was tested and optimized with human and animal macrophage growth in vitro. It is formulated for use with 5% CO₂ and 95% air in a humidified incubator. The medium contains essential and non-essential amino acids, vitamins, organic and inorganic compounds, hormones, growth factors, trace minerals), supplemented with macrophage growth supplement, antibiotics, and fetal bovine serum.

Storage Condition

Store the basal medium at -20°C until use. The complete cell culture medium can be kept at 4°C for 6 weeks. Protect from light.

To assure sterility, the prepared medium should be filtered with a 0.2 um filter after 2 weeks or if there is concern that sterility was compromised during the supplementation process.

Shipping

Dry ice.

Authorized Uses of *Cell Biologics* Products

Macrophage Media from *Cell Biologics* are distributed for internal research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' cells or products from the purchaser to other markets, organizations or individuals is prohibited by *Cell Biologics*, without the company's written consent. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.