

## **Trypsin-EDTA 0.25% Solution (1X) - 100ML**

Catalog Number: 6914

### **Description**

Trypsin solutions are used widely for dissociation of tissues and cell monolayers and we can deliver solutions that will perform consistently, batch after batch.

The biochemical assays performed on Trypsin 1:250 determine both trypsin specific-activity at the level of certain co-purified enzymes that influence cell removal and viability. All trypsin solutions are prepared from porcine-parvovirus-tested and mycoplasma-tested materials.

### Formulation

Contains 0.25% trypsin and 1mM EDTA in Hank's Balanced Salt Solution without calcium and magnesium.

Sterile filtered.

### **Storage Condition**

Store the solution at -20°C.

### **Shipping**

Ice pack.

### **Authorized Uses of *Cell Biologics* Products**

Trypsin-EDTA 0.25% Solution from *Cell Biologics* are distributed for internal *in vitro* research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' cells or products from the purchaser to other markets, organizations or individuals is prohibited by *Cell Biologics*, without the company's written consent. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.