

**Collagen Type I, Rat Tail Solution - ready to us**

Catalog Number: 6953

**Description**

Collagen Type I, Rat Tail Coating Solution is designed to coat cell culture flasks or dishes. Coat sterile culture dishes or flasks with Collagen Type I, Rat Tail Coating Solution for 2 hours, then aspirate the excess solution and wash dishes with PBS once before seeding cells. The use of coated-dishes improves cell attachment for certain types of primary cells.

**Storage Condition**

Store the solution at -20°C.

**Shipping**

Ice pack.

**Authorized Uses of *Cell Biologics* Products**

Collagen Type I, Rat Tail Solution from *Cell Biologics* are distributed for internal *in vitro* research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' cells or products from the purchaser to other markets, organizations or individuals is prohibited by *Cell Biologics*, without the company's written consent. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.