

**Endothelial Cell Basal Medium (Without Glucose and Phenol Red) – 500 ML**

Catalog No: GPF1168b

**Description**

Endothelial Cell Basal Medium without glucose and phenol red is a medium designed for the culture of endothelial cells. It was tested and optimized with human or animal endothelial cell growth and proliferation *in vitro*. It is formulated for use with 5% CO<sub>2</sub> and 95% air in a humidified incubator.

The basal medium contains essential and non-essential amino acids, vitamins, organic and inorganic compounds, hormones, growth factors, trace minerals.

**Storage Condition**

Store the basal medium at 2-8°C for 8 months. Protect from light.

**Shipping**

Ice pack.

**Authorized Uses of *Cell Biologics* Products**

Cell Culture Medium from *Cell Biologics* are distributed for internal research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' Cells or Products from the purchaser to other markets, organizations, or individuals is prohibited by *Cell Biologics*. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.