

**RPMI-1640 Media; With L-glutamine and HEPES Buffer**

Catalog No. M6317

**Description**

RPMI-1640 was developed in 1966 by Moore and his co-workers at Roswell Park Memorial Institute (hence the acronym RPMI). It is based on the RPMI-1630 series of media utilizing a bicarbonate buffering system and alterations in the amounts of amino acids and vitamins. While it was originally formulated to support lymphoblastoid suspension culture, it has been proven to support a wide variety of anchorage-dependent cells. It has a wide range of uses, including the culture of fresh human lymphocytes, fusion protocols, and in the growth of hybrid cells.

**Storage Condition**

Store the basal medium at 2-8°C. Store cell growth supplement, fetal bovine serum (FBS) and antibiotics at -20°C. The complete cell culture medium with Supplement Kit can be kept in 4°C for two months. Protect from light.

To assure sterility, the prepared medium should be filtered with a 0.2 um filter after 2 weeks or if there is concern that sterility was compromised during the supplementation process.

**Shipping**

Ice pack.

**Authorized Uses of *Cell Biologics* Products**

Cell Medium from *Cell Biologics* are distributed for internal *in vitro* research purposes only. Our products are not authorized for human use, for *in vitro* diagnostic procedures, or for therapeutic procedures. Transfer or resale of any *Cell Biologics*' cells or products from the purchaser to other markets, organizations or individuals is prohibited by *Cell Biologics*, without the company's written consent. *Cell Biologics*' Terms and Conditions must be accepted before submitting an order.