

BALB/c Mouse Primary Splenocytes

Catalog No. BALB-6321

Suggested Medium

M5567 Complete Cell Culture Medium/w Kit - 500 ML

Product Description

BALB/c Mouse Primary Splenocytes are derived from the spleen of pathogen-free laboratory mice. Each vial contains 10×10^6 cells and is delivered frozen. Cells are negative for bacteria, yeast, fungi, and mycoplasma. Cells can be used for designed experiments under the cell culture conditions specified by Cell Biologics. These cells do not proliferate in culture and cannot be passaged.

Storage

Cryopreserved cells will be shipped overnight on dry ice. Upon arrival, please immediately transfer the frozen cells to liquid nitrogen (-180°C) until ready for use. Suspension cells can be shipped in 50 ml conical tubes upon request. Primary cells should never be stored in a -20°C or -80°C freezer.

Authorized Uses of Cell Biologics' Products

BALB/c Mouse Primary Splenocytes from Cell Biologics are distributed for internal research purposes only. Our products are not authorized for human use, for in vitro diagnostic procedures or for therapeutic procedures. Transfer or resale of any Cell Biologics' cells or products from the purchaser to other markets, organizations or individuals is prohibited by Cell Biologics, without the company's written consent. Cell Biologics' Terms and Conditions must be accepted before submitting an order.

Disclaimer

Investigators should handle the cells that they receive from Cell Biologics with caution and treat all Cells as potential pathogens, since no test procedure can completely guarantee the absence of infectious agents.

Warranty and Liability

Cell Biologics' guarantee applies only to your purchase of Cell Biologics' cells with Cell Biologics' Media and Coating Solution for appropriate cell culture and cell testing following Cell Biologics' online protocols within 35 days from the date of product delivery.